## Wireless Body Composition Smart Scale User Manual

The Wireless Body Composition Smart Scale is body analyzer that measures weight, and it uses bioelectrical impedance analysis (BIA) technology to estimate body fat percentage and other body composition data. It is intended for home use only.

### Specifications

Dimensions: 260 x 260 x 24 mm (10.2 x 10.2 x 0.9 in)  
Product weight: 1.0 kg (2.2 lb)  
Batteries: 2 x 1.5V AAA  
Measurement capacity: 0.2 ~ 180 kg (0.5~400 lb)

### Display Description

	(no display)	
Normal working condition	No battery / battery depleted	Wireless not connected
Low on power	Overload	Calibrating

### WARNING

Please read carefully before use

- Do not use if you have a pacemaker, or artificial lung, or wear an ECG, or any kind of internal medical device.
- Do not use on slippery floor; do not use when your body is wet.

### IMPORTANT NOTES

- Your body composition data measurements may not be accurate when using the scale immediately after heavy physical activities.
- The scale is best suited for adults over age of 18. (Children between 0-9 may avail themselves of weight and BMI data, teenagers 10-17 of weight, BMI and bodyfat percentage, but other data are best served as reference).
- The Smart Scale should not be used to diagnose or treat any medical condition. Always consult your physician or medical professionals to determine what is ideal for you.

### Power On Scale

Your Smart Scale is shipped with 2 AAA batteries. Install the batteries inside the battery compartment and your scale is power on.

With batteries installed, you can check your body weight by stepping on it. However, to measure other body composition data such as bodyfat percentage, you need to use your Smart Scale together with your smartphone. To learn how to install mobile application program on your mobile devices, see next section.

### APP Install and Measurement


The Smart Scale connects wirelessly to your mobile phones. To do so, your mobile devices need to install our recommended mobile application program. The easiest way to download the mobile APP is to use your phone to scan the QR code on following page. Follow the instructions to experience the complete feature of the scale.

### Measurement Tips

Read these tips to help you get the most accurate measurements:


- This product employs BIA (bioelectrical impedance analysis) technology. To have accurate and consistent measurement result, please use it at the same time during the day and under same conditions as much as possible. Please avoid taking measurements:
  - During or immediately after strenuous exercises
  - After a heavy meal
  - Immediately after a shower, bath or sauna
  - With heavy alcohol consumption
  - During pregnancy
  - During fever
- During the test, please be barefooted and stand straight. Wearing socks may not get you the accurate result. Dirty feet may also affect accuracy of the measurement.
- Weight unit can be changed from pound (lb) to kilogram (kg), or vice versa, via the mobile APP program. The height unit will also change accordingly (meter vs. foot/inch).

- Measure weight correctly: place the scale flat on hard floor and reset it before taking measurement. The weight should display 0.0 kg/lb. before stepping on the scale.

Use it on hard surface

- Measure bodyfat and other body composition parameters correctly: Stand straight on scale, and have both feet uniformly touched the electrodes as figure below.

Conductive areas

### Step 1: Download the APP

Scan the QR Code or search APP at the following link: [app.cubitt.co/scale](http://app.cubitt.co/scale) to download mobile APP


Register and fill in your personal data


Install battery cells to the scale


### Step 2: Pairing scale with the APP

Turn on bluetooth of phone Enter my device/Add new device/Pair device

Step on scale


Pairing scale to APP

The pairing is automatic once you step on the scale, otherwise you should connect the scale manually in the **Setting menu of the APP**.

**Note:** Device running Android 6.0 will require enabling Location Access for the pairing to work.


### Step 3: Measure and check the body analysis report

Quick-step on the scale for calibration (can be done anytime)

Step on again until the display number is stabilised

Check the body analysis report in the APP

Number flashes twice indicates weight test finishes. Number 0 rolls to the right then stabilises indicates the bodyfat test has also completed.

### What's Included

Your Smart Scale Package includes:

- Smart Scale
- Battery: 2 x 1.5V AAA
- User Manual

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